

Induced pluripotent stem cell (iPSC)-derived hepatocytes demonstrate functional Phase I and Phase II drug metabolism and can accurately predict drug-induced liver injury (DILI)

Malek Haddad¹, Magdalena Lukasiak¹, Gemma Gatti¹, George Kiloh¹, Samuel Chung¹, Ka Chung¹, Chloe Robinson¹, Bomi Jung², Bhavik Chouhan², Dominic Williams³, Lia Panman¹, Carlos Gil¹, Nikolaos Nikolaou¹

¹DefiniGEN Ltd., Cambridge, United Kingdom; ²AstraZeneca, Gothenburg, Sweden; ³AstraZeneca, Cambridge, United Kingdom

B

INTRODUCTION

Liver disease is a rising cause of mortality worldwide. Despite careful approaches taken during drug development, most therapeutics still fail to reach clinical trials due to the lack of translatability between preclinical models and the clinic. Primary human hepatocytes (PHH) and hepatocellular cancer cells are currently the predominant pre-clinical liver models in the industrial workflow. However, they come with limitations, including short life-span, rapid loss of function, and malignant origin. We have hypothesized that induced pluripotent stem cell (iPSC)-derived hepatocytes can overcome these limitations, offering a sustainable in vitro platform that is more relevant to human liver

Ulti-HEP SHOW COMPARABLE PHASE | AND PHASE || EXPRESSION OF DRUG METABOLISM ENZYMES TO PRIMARY HUMAN HEPATOCYTES



Figure 3: A) mRNA expression levels of Phase I CYP450 genes in hepatocellular carcinoma HepG2 cells, DefiniGEN Ulti-HEP, and primary human hepatocytes (PHH). mRNA data were normalized to the housekeeping gene 18S rRNA and are presented as mean±SEM of 3-4 independent experiments. B) mRNA expression levels of UDP-glucoronosyltransferases in Ulti-HEP and PHH. mRNA data are presented as raw counts following bulk RNA-sequencing analysis. For PHH data, cells from 1-3 independent donors were used.

physiology, amenable, and expandable for large-scale drug screening and predictive toxicology studies.

MATERIALS AND METHODS

Healthy iPSCs were differentiated towards hepatocytes (Ulti-HEP) in 2D format for six weeks. Following iPSC differentiation, characterization of Ulti-HEP functionality was assessed, including liver maturity marker expression by immunocytochemistry, urea synthesis and secretion by western blotting and biochemical assays, and Phase I and Phase II drug metabolism-associated gene expression by qPCR, alongside that of hepatocellular carcinoma HepG2 cells and PHH. Suitability of Ulti-HEP to predict drug-induced liver injury (DILI) was finally evaluated by cell viability assays (ATP endpoint) following a 48 hour-incubation with drugs with known DILI liability, alongside metabolite analysis of Ulti-HEP treated with acetaminophen and diclophenac and analysed by Liquid Chromatography/ Mass spectrometry analysis.

Ulti-HEP DEMONSTRATE COMPARABLE HEPATOCYTE MATURITY MARKER EXPRESSION TO PRIMARY HUMAN HEPATOCYTES

Ulti-HEP SHOW DEMONSTRATE PHASE | AND PHASE II DRUG METABOLISM

Figure 1: Representative immunocytochemistry pictures showing protein expression of the hepatocyte maturity markers albumin (red), alpha-1-antitrypsin (A1AT; green), HNF4A (green), and AFP (red) in liver carcinoma HepG2 cells, DefiniGEN Ulti-HEP, and primary human hepatocytes (PHH). Nuclei were counterstained with DAPI. Data are presented as mean±SEM of n=3-4 independent experiments.

Figure 4: A, C) Simplified schematics demonstrating acetaminophen and diclofenac metabolism within human liver. **B**, **D**, **E**) Wild-type Ulti-HEP media metabolite analysis following treatment with acetaminophen or diclofenac, as measured by Liquid-Chromatography/ Mass spectrometry analysis. Ulti-HEP in (D) were treated with increasing diclofenac concentrations and compared against a 3D PHH model treated with 10 µM diclofenac. Data are presented as mean±SD of n=3 technical replicates.

Ulti-HEP DEMONSTRATE A LONG-TERM FUNCTIONAL UREA CYCLE

Figure 5: A) Cell viability (ATP content) in DefiniGEN Ulti-HEP following 48 hours of treatment with increasing concentrations (0-250 μ M) of 7 compounds with known DILI liability. Cell viability data were normalized to positive (250 µM chlorpromazine) and negative controls (0.2% DMSO) and are presented as mean±SD of n=3 technical replicates. B) The ratio of cellular ATP IC50 to total plasma Cmax (RA/C) was calculated and plotted; the RA/C values of the "high-DILI-concern" drugs cluster below 100, whereas the RA/C values of the "no-DILI-concern" drugs cluster above 100. The RA/C values of the "medium-DILI-concern" drugs cluster between the two other groups.

DISCUSSION

Figure 2: A) The urea cycle which metabolises ammonia to urea, highlighting the enzymes involved in green boxes. B) Protein expression levels of the urea cycle enzymes OTC, ASS1, ASL, CPS1, and ARG1 in liver carcinoma HepG2 cells and DefiniGEN Ulti-HEP. C) Urea secretion in DefiniGEN Ulti-HEP over a period of three weeks. Data are presented as mean±SEM of n=3-4 independent experiments.

To conclude, we have shown that Ulti-HEP offer a superior in vitro liver model compared to HepG2 cells, with greater liver functionality that is comparable to primary human hepatocytes. These data alongside the expansion capacity and amenability of Ulti-HEP using CRISPR-based gene editing showcase the spectrum of opportunities iPSC-derived hepatocytes can offer in the fields of disease modelling, large-scale drug screening, genotoxicity, and hepatic safety.